

Article Datasheet

Art.No.: 81792

Refill Bottle for Eye Wash Station

WITH STERILE SODIUM CHLORIDE SOLUTION (0.9%)

Article Description:

- Bottle with sodium chloride solution (0.9%) for use to protect eye after incident
- suited with ergonomically formed eye- and dust-protection
- 200ml bottle for mobile use, like first-aid-kit and cars
- 500ml bottle for one-time-use or as refill for eye wash station
- time of use: about 5 minutes
- Packing Unit: 200ml (.../200) = 10, 500ml (.../500) = 12



Catalogue Page: T 115